

Aerosol Particle Density Determination Using Light Scattering in Conjunction with Mass Spectrometry

Eben S. Cross¹, Paul Davidovits¹, John T. Jayne², David K. Lewis^{2,3}, Manjula Canagranata², James Allan⁴, Doug Worsnop², and Tim Onasch²

¹Boston College, Chestnut Hill, MA 02467

²Aerodyne Research Inc. Billerica, MA 01810

³Connecticut College, New London, CT 06320

⁴The University of Manchester, Manchester M60 1QD UK

ASP Meeting

Boulder, CO


October 25 – 27, 2006

Acknowledgements: NSF grant No. ATM-0525355. DOE grants No. DE-FG02-05ER63995 and No. DE-FG02-05ER84268. ESC funded by NASA Earth System Science graduate student fellowship


Features of the LS-AMS

- All particles that impact the vaporizer surface ($d_p > 180$ nm), first pass through the laser beam and scatter light – providing per particle measurements of scattered light (R_{LS}), vacuum aerodynamic diameter (d_{va}) and chemical ion signal (R_{MS}). (Panel 1 and 2)
- Correlated R_{LS} (optical) and R_{MS} (chemical) measurements provide an in situ measurement of the collection efficiency of the AMS. (Panel 2)
- $R_{LS} - d_p$ calibration curve allows the determination of an optical diameter (d_o). (Panel 3)
 - Using a single species calibration curve for spherical particles d_p is determined to +/- 10% accuracy across the refractive index range of $1.41 < n < 1.60$. (Panel 4)
- Combination of d_o and d_{va} provides a single particle density determination, $\rho_{LS} = (d_{va}/d_o)$.
- The LS-AMS per particle density determination allows the mixing state of the ambient aerosol particles to be analyzed.
- Single particle density (ρ_{LS}) distributions obtained for ambient aerosol particles at Chebogue Point, Nova Scotia are highly correlated with the chemical composition-based density (ρ_{cc}) and show that the ambient particles are internal mixed. (Panel 6 and 7)


Schematic of the LS-AMS




Correlated R_{LS} and R_{MS} Signals



Experimental vs. Theoretical Scattering Response of the LS-AMS System



Effect of Refractive Index on R_{LS}




*Note: n = real component to complex refractive index. Absorption component for oil droplets ($k \sim 0$)

LS-AMS Deployment: Ambient Aerosol Analysis



July 11 – August 15, 2004
Chebogue Point, Nova Scotia




- Residence time of the aerosol particles in the atmosphere ~ days

6

Real-Time Determination of Chemical Composition,

$$\rho_{cc} \text{ and } \rho_{LS}$$



Correlation Between the two Density Measurements:

$$\rho_{cc} \text{ VS } \rho_{LS}$$

