

Direct Aerosol Absorption Measurements using Photothermal Interferometry

Arthur J. Sedlacek III and Jeonghoon Lee
Brookhaven National Laboratory
Atmospheric Sciences Division
sedlacek@bnl.gov

Introduction

One of the goals of the atmospheric community is to better quantify the role of aerosol absorption on the radiative balance of the atmosphere. Black carbon (BC) is important since it absorbs light thereby contributing to global warming. However, the level of uncertainty regarding the contribution of aerosols (BC & BC-containing particles) to radiative forcing is poorly understood.

The goal of this project is to develop a method for directly measuring aerosol absorption that is capable of aircraft deployment.

Photothermal Methods for Measuring Aerosol Absorption

Photothermal Spectroscopy (PTS), which relies on the thermal dissipation of spectrally absorbed energy for its signal, offers the highly desirable trait of directly measuring aerosol absorption

Two PTS-based approaches are currently under development or in use for aerosol absorption:

Photoacoustic Spectroscopy (Arnott & co-workers, 1998, 1999, 2005; Lack, et al., 2006)

&

Photothermal Interferometry (Sedlacek, 2006)

The hallmark of any photothermal-based technique is that it directly measures aerosol absorption without interference from aerosol scattering

Photothermal Interferometry

Interferometry is the study of the interference patterns created by combining two sets of coherent light waves that have a constant phase relation and same polarization. This very sensitive technique can be used to follow the thermal dissipation that characterizes PTS.

$$E_r = A_r e^{i(\omega_r + \phi_r)} \quad E_p = A_p e^{i(\omega_p + \phi_p + \Delta\phi(t))}$$

$$\text{where } \Delta\phi(t) = m \sin(\omega_{\text{mod}} t)$$

$$\text{Intensity} = (E_r + E_p)(E_r + E_p)^* \dots$$

$$I(t) = [A_r^2 + A_p^2 + 2A_p A_r [\cos(\phi_r - \phi_p) \cos \Delta\phi(t) - \sin(\phi_r - \phi_p) \sin \Delta\phi(t)]]$$

$$\text{In the limit of } \Delta\phi \ll 1: \phi_r - \phi_p = \pi/2$$

$$I(t) \propto C \Delta\phi(t)$$

Photothermal Interferometry: Performance Estimate

$$\begin{aligned} \Delta\phi_{\text{min}} &= 1.6 \times 10^{-8} \text{ rad} \\ r &= 0.5 \text{ mm} \\ \rho &= 1.205 \text{ kg m}^{-3} \\ C_p &= 1005.4 \text{ J kg}^{-1} \text{ K}^{-1} \\ \omega_{\text{mod}} &= 100 \text{ Hz} \\ \lambda &= 632.8 \text{ nm} \\ l &= 5 \text{ cm} \\ n &= 1.000292 \\ P_{\text{exc}} &= 250 \text{ mW} \\ T_o &= 293 \text{ K} \end{aligned}$$

$$\alpha_{\text{min}} = \frac{\Delta\phi_{\text{min}} 2r^2 \rho C_p \omega_{\text{mod}} \lambda T_o}{\pi(n-1)P_{\text{exc}}}$$


$$\alpha_{\text{min}} \sim 5 \times 10^{-8} \text{ m}^{-1} = 0.05 \text{ Mm}^{-1} = 50 \text{ Gm}^{-1}$$

Relevant PTI Studies


Lin and Campillo (1985) used PTI to measure $(\text{NH}_4)_2\text{SO}_4$ absorption using a Mach-Zehnder interferometer

Detection sensitivity of 50 Mm^{-1} ($\lambda_{\text{abs}} = 1086 \text{ cm}^{-1}$) – sensitivity was limited by bgkd water vapor absorption at this wavelength

More recently Davis and co-workers (1999) using a modified Jamin interferometer reported a sensitivity of 37 Gm^{-1} for NH_3 ($\lambda_{\text{abs}} = 1085 \text{ cm}^{-1}$).




Folded Jamin Interferometer Design offers superior optical stability



- Insensitive to rotation and tilt
- Simple positional alignment
- Common path enables common mode rejection providing increased robustness, stability (mechanical vibrations) and sensitivity
- Compact size
- Maximal pathlength for given size (inherent double-pass configuration)


AIMS folded Jamin interferometer manufactured by Interferomet Inc.
<http://www.interferomet.com>

BNL Folded Jamin-based PTI Design




Glancing configuration eliminates all spurious signals

PTI unit is designed to fit into standard instrument rack




Subsystem Performance Characteristics



Servo-control requirements: Maintain system at quadrature by correcting for slowing varying drifts (e.g., thermal) but passing high-frequency signal onto LIA.

Drifts on the order of <100 Hz are corrected by feedback circuit





Use Lissajous circle to optimize unit

Unit Calibration Measurements: NO₂ @ 514 nm

Calibration standard (410 ~ 600 nm)

Relevant PSD parameters:
Detection bandwidth: 7.8 mHz
Pathlength: 5 cm
RMS noise floor: 300 nV



Measured absorption cross-section compares favorably with previous literature values:


$$\sigma_{\text{abs}}^{\text{PTI}} = 2.2 \pm 0.08 \times 10^{-18} \text{ cm}^2 \text{ molec}^{-1}$$

$$\sigma_{\text{abs}}^{\text{lit}} = 1.9 \pm 0.18 \times 10^{-18} \text{ cm}^2 \text{ molec}^{-1} \text{ (Average of 4 different studies*)}$$


*Vandaele et al., JAC, 25 1996; Davidson et al., JGR, 93, 1988; Harder et al., JGR, 103, 1997; Vandaele et al., JQSRT, 59, 1998

Evaluation of PTI unit Towards Measurement of Aerosol Absorption: Nigrosin


Nigrosin is a polyaniline-based black pigment. Its broad, featureless absorption spectrum in the visible makes it an ideal BC surrogate for instrument evaluation.




Size-selected nigrosin aerosols were generated using TSI constant output atomizer and DMA and counted using CPC.



Using our calibration standard (NO₂), the absorption coefficient (Mm⁻¹) of nigrosin aerosols as a function of particle diameter can be determined. Representative data for three different particle diameters is shown below.




Measurement of the PTI signal as a function of nigrosin number density gives a slope that is equal to the absorption cross-section (cm²/particle). Shown below are some measurements of σ_{abs} for various nigrosin particle diameters for different manufacturers superimposed on the recently published results by Lack et al. (2006).




Measurements on non-absorbing aerosols: (NH₄)₂SO₄

Despite the concentration being 4-5 orders-of-magnitude higher than that used in the nigrosin experiments, no signal was observed, demonstrating that PTI does not suffer from scattering-related interference effects



Measurement of Ambient Aerosols

Experiments have recently been initiated to measure ambient aerosols at BNL using PTI (B_{abs}) and an integrating nephelometer (B_{scat}). The ability of the PTI technique to directly measure B_{abs} enable accurate measurement of SSA values, especially very near unity.



Publications

Sedlacek, A. J., "Real-time detection of ambient aerosols using photothermal interferometry: Folded Jamin interferometer" Rev. Sci. Instr. 77 064903 (2006)

Acknowledgements

Authors would like to thank Dr. J. Wang for use of his DMA and CPC.