


Collaborative Aerosol Research

Boston College – Aerodyne

Onasch^{1,2}, Davidovits¹, Cross¹, Hings¹, Wrobel¹, Slowik³,
Han, Trimborn², Jayne², Williams², Kroll², Canagaratna²,
Kolb², Worsnop^{1,2}


¹Boston College; ²Aerodyne Research, Inc.;
³ University of Toronto



**DOE Atmosphere Science Program
FY 2008 Science Team Meeting
Annapolis, MD
February 25-27, 2008**





Aerosols in the Atmosphere




FOCUS: Characterize chemical and physical properties and atmospheric transformations of aerosol particles through laboratory studies and extractive sampling techniques

- **Boston College Laboratory: Carbonaceous aerosol generation, modification, and characterization**
 - Sooting flame generation
 - DMA-AMS, SMPS, and CCN chemical and physical characterization
 - Black carbon instruments comparison studies (e.g. MAAP, SP2, PAS)
- **Field studies**
 - Mexico City (2002, 2003, 2006)
 - NEAQS (2004)
- **Instrument development and testing**
 - LS-QAMS
 - LS-C-TOF-AMS
 - SP2-AMS




Laboratory: Flame Apparatus



- Premixed flame: fuel-to-air (equivalence) ratio from <1 to >5 (sooting flame above 1.8)
- Controlled sampling system with reproducible particle generation
- Polydisperse or DMA size-selected soot particles
- Sample flow rates of 0.2 – 4 lpm provided to multiple instruments
- Thermal denuder and Condenser to modify generated soot particles

Monodisperse Soot Generation




Operating Parameters

$$1.8 < \phi < 5.0$$



$$30 < D_{\text{mobility}}(\text{nm}) < 500$$

$$15 < D_{\text{pp}}(\text{nm}) < 55$$

$$0 < \text{Coatings}(\text{nm}) < 150$$



Soot Characterization, Coating, and Denuding Experiments



- Particle shape affects size measurements and derived mass measurements
- Results from premixed laboratory flame experiments show behavior similar to diesel soot particles

Laboratory and Field: Chemical flow tube


Introducing the concept of Potential Aerosol Mass (PAM)

E. Kang¹, M. J. Root¹, D. W. Toohey², and W. H. Brune¹

¹Department of Meteorology, Pennsylvania State University, University Park, PA 16802, USA


²Atmospheric and Oceanic Sciences, University of Colorado, CO 80309-0311, USA

Atmos. Chem. Phys., 7, 5727–5744, 2007



- Collaboration with Prof. Brune in PAM characterization and employment
- Characterization of Highly Oxidized Secondary Organic (SOA) Coatings
- Specifically looking at AMS mass spectra of organic compounds as a function of oxidation (e.g. OOA2 to OOA1 progression)




CCN (cloud activation) experiments



- Particles composed of soluble (e.g. ammonium sulfate), slightly soluble (e.g. adipic acid), and insoluble (flame soot) materials

*****poster*****

Adipic Acid, Ammonium Sulfate, Soot



- Measuring CCN over large size and SS ranges
- Particle phase is very important for slightly soluble compounds
- Köhler theory works well mixed adipic acid and sulfate particles
- Insoluble soot core acts as support structure for small amounts of adipic acid coating, but larger soot particles may not get fully coated


****poster****

Instrument Development: LS-C-ToF-AMS

MS Mode


PToF Mode

LS-SP Mode




- Practical single particle aerosol mass spectrometer
- Alternating with measurements of ensemble chemistry and size

Mexico City Field Deployment



<http://moe2.org/megaioites/fieldcampaign2006/>





- LS-C-TOF-AMS operated for 74 hours at the T1 site NE of Mexico City
- Measured single particle mixing states while contemporaneously measuring ensemble chemistry and size

in collaboration with Liz Alexander (EMSL, PNNL)

****poster****


Direct Observations of Condensation onto Fresh Soot Particles



- Single particle detection of secondary oxidation products (NO₃ and SOA) condensing on primary HOA particles

*****poster*****

Instrument Development: SP2-AMS




- Absorbing particles (coating and core) vaporize in laser
- Vapor is ionized by electron impact ionization
- Detection of the ions by Time-of-Flight mass spectrometry

in collaboration with Greg Kok at DMT


Chemistry, Mass, and Size of Coatings and Core

Oil Lamp Soot Particles



- Obtain chemical information on elemental carbon clusters and compounds coating the soot cores
- Mass signals are linear with respect to coating thickness and core size
- PTOF show size and particle mass signals increasing with particle coatings


225 nm Glassy Carbon Spheres uncoated and coated with DEHS oil




****poster****

Inter-Comparison of Instruments Measuring Black Carbon Content of Soot Particles


Photoacoustic Spectrometer (PAS)




Multi-Angle Absorption Photometer (MAAP)




Single Particle Soot Photometer (SP2)



Aerosol Mass Spectrometer (AMS)



Instrument Comparisons



- Carbon spheres, uncoated flame soot, coated (liquid and solid) flame soot
- Flame equivalence ratios of 1.8 to 2.3
- Particle shape and solid coatings affected optical measurements

Planned Black Carbon Measurements



- Summer 2005 – initial BC study
- Summer 2008 – second Boston College laboratory study
 - Improved Flame Apparatus
 - Soot with equivalence ratios 1.8 to 5.0 (varying D_{pp})
 - Thermal denuder
 - Condenser for organics and inorganics (sulfuric acid)
 - RH control in sampling lines
 - Instruments: SP2-AMS (Aerodyne/DMT), Photothermal Interferometric Absorption Spectrometer (Sedlacek/BNL), MAAP (Aerodyne), Photo-Acoustic Soot Spectrometer (DRI/LANL), SP2 (NOAA/DMT), OC/EC (?),

Summary

- Laboratory studies generating, modifying and characterizing carbonaceous soot particles
- CCN of mixed particles composed of soluble, slightly soluble, and insoluble components
- LS-C-TOF-AMS instrument development as a single particle aerosol mass spectrometer
- SP2-AMS instrument development for carbon particle chemistry, mass, and size measurements
- Black carbon instrument inter-comparison studies


SP2-AMS Sampling

Vaporization and Ionization Mechanism



- Coatings evaporate first at relatively low temperatures (<600oC) potentially dependent upon vapor pressures
- Core evaporates last at high temperature (>1000oC) while incandescing (SP2)
- 5-20 microsecond evaporation time

Aerosols in the Atmosphere



C.E. Kolb, *Nature*, 597-598, 6 Jun 2002