



Cloud Microstructure via the Fast Time-Resolved Aerosol Collector Fast TRAC.....

...with Holographic Imaging

**Xiao-Ying Yu, Ali Hashim,
Martin Iedema, and James Cowin**

**Atmospheric Sciences, Chemical Sciences
Pacific Northwest National Laboratory
Richland, WA**

Research is supported by NOAA & DOE.

Cloud Microstructures ≤ 1 m


Want to know the aerosols at this resolution





Aircraft flies at 150 m/s

Need time resolution 1 m/150 m/s
= 6 ms (!!!!!)


What is TRAC? - Time-Resolved Aerosol Collector



- Uses an impactor
- ~ 600 TEM samples
- Flow rate: 1 l/min
- Time resolution: ≥ 1 min*
- Applications: Off-line analysis:
 - particle hygroscopicity, morphology, composition..




Fast TRAC Solution




- Observe the particles DURING their collection with video microscope
- See ≥ 100 nm particles
- ~ 4 ms time resolution


Results - movies of 200 nm lab particles




Xiao-Ying Yu, [James Cowin PNNL](#)



Difference
photo is black
except for
diffraction-limited
image of 200 nm
particle




Particle Locations to ± 0.1 micron




Real-time Optical Sizing !!!!

Xiao-Ying Yu, [James Cowin PNNL](#)

Optical and STEM Photos




Optical particle map and times




SEM/EDX map and elemental analysis

Fast Framing



Frame 153




Frame 154

New 100 nm
Particle Arrives

- 240 frames/sec
- 4 ms per frame
- See ≥ 100 nm particles arriving

Fast TRAC Features

Size of new Fast-TRAC



- Real-time particle optical sizing
- 4 ms time resolution
- Extensive off-line analyses
- **Good for cloud microstructures**
- **And plumes**


We need longer depth of field....

Substrates not always perfectly flat

Severe vibrations may shift , flex film





To my surprise, imaging works well in brightfield illumination too, even for 100 nm particles...

1) Improve depth of field by incoherent well collimated brightfield illumination



2) Improve depth of field by post processing for coherent holographic brightfield illumination

...as used in some snowflake imagers, undersea imagers



Future Work

- Make it field-portable
- Deploy Fast TRAC in field campaigns
- Collaborations

PNNL's Fast TRAC for cloud microstructures and plumes